


CARD–FISH and Microautoradiography Protocol for Bacteria and Archaea



by Gerhard Herndl Lab @ www.microbial-oceanography.eu

No liability shall be taken for direct, indirect, consequential or incidental damages arising from the use, results of use, or inability to use this protocol. The protocol has been thoroughly checked, however, if you have suggestions for improvements please let us know (thomas.reinthal@univie.ac.at).

Sample Fixation

1. Fix samples with paraformaldehyde or 35% formaldehyde to final concentration of 2-4% and fix for at least 1 but not longer than 18h.
2. Filtrate sample (for $\sim 10^6$ cells/ml 10ml will do) on white polycarbonate (25 mm 0.2 μ m) filter with cellulose nitrate support filter (0.45 μ m).
3. After sample filtration, wash with 5-10 ml of MQ
4. Air-dry filters
5. Store at -20°C until processing. Filters can be stored frozen for several months.

Embedding

1. Warm up 0.1% Agarose and pipette 30 μ l on petri dish. Dip filter with both sides in Agarose and place filters shiny side down on drop
2. Let filters dry at 37°C in hybr. oven without lid for 10-15min.
3. Pour ethanol (95%) in petri dish and carefully remove filters
4. Dry filters, now they can be stored at -20°C

Permeabilization

1. Prepare 10ml of permeabilization mix with **either** lysozyme (Eub) **or** proteinase K (Arch); after preparation put on ice

Stock reagent	Volume (μ l)	Final
Lysozyme	100mg	10mg/ml
Proteinase K	2check	10.9mg/ml
1M Tris-HCl	1000	0.1M
0.5M EDTA	1000	0.05M
MQ water	8000	

2. Pour lysozyme (proteinaseK) in petri dish and place filters upside down into it, incubate for 1h at 37°C in hybr. oven
3. After 1h wash filters in MQ (1x for Lysozyme, 3x for Proteinase K)
4. Place filters for 20-25min. in 0.01M HCL at RT
5. Wash filters 2x in excess MQ
6. Shortly dip them into Ethanol (95%)
7. Dry filters, now they can be stored at -20°C .

Hybridization

1. Prepare hybridization buffer:
2. Cut 25mm filters in up to 8 sections and place them into 0.7ml Eppi (labeled).
3. Mix 300 μ l hybridization buffer + 15 μ l HRP probe. Probe to buffer ratio: 1:20. Freeze probes just once, after thawing store them at 4°C for up to 0.5 years.
4. Hybridize at 35°C for 12h – 15h (low activity), stick Eppis around rotor, darken them and shake slowly
5. Prepare washing buffer:

Washing Buffer for 35°C and 55% Formamide in Hybridization buffer:

Stock reagent	Volume (μ l)	Final
5M NaCl	30*	13mM
1M Tris-HCl	1000	20mM
0.5M EDTA	500	5mM
MQ water	48420	
10% SDS	50	0.01%

Washing Buffer for 35°C and 20% Formamide in Hybridization buffer:

Stock reagent	Volume (μ l)	Final
5M NaCl Arch	1350*	145mM
1M Tris-HCl	1000	20mM
0.5M EDTA	500	5mM
MQ water	47100	
10% SDS	50	0.01%

1. Pre warm buffer at 37°C.
2. After hybridization quickly transfer filter sections 15 min in washing buffer, afterwards pour in to buechner funnel
3. *chemical corrections applied (mM -100 μ L)

Amplification

1. Prepare PBS-T-Mix:

Stock reagent	Volume (μ l)	Final
1xPBS	49750	
100% Triton X100	25	0.05%

PBS should have pH of 7.4 -7.6 (better peroxidase turnover)

1. Pick filters and incubate in PBS-T-Mix at RT for 10-15min.
2. Prepare substrate mix with prepared amplification buffer (AMP) and 30% H2O2
3. Dilute 30% H2O2 to final concentration of 0.0015% (Step A to B; all Volumes in μ l)

	AMP	H2O2	Tyr	Ratio	
A	200	1 from 30%			
B	493	5 from A	2.5	1:200	A555
B	493	5 from A	3.4	1:150	A488
B	493	5 from A	5.0	1:100	

1. Dab filters on blotting paper, place sections on top of each other and place in 1.5ml Eppi
2. Incubate in substrate mix for 10min up to 45min. at 37°C in hybr. oven in the dark
3. After incubation in Substrate mix dab filters on blotting paper
4. Wash in 50ml PBS-T mix at RT for 10min in the dark
5. Wash in MQ
6. Wash in Ethanol 95%
7. Dry and mount in DAPI mix

Hybridization Buffers - 10ml:

Store at -20°C for up to 1 year; if in use keep on ice !

Buffer with 55% Formamide:

for EUB, ALF, BET, GAM, CF

Stock reagent	Volume (µl)	Final conc.
Dextran Sulfate	1g	10%
5M NaCl	1800	900mM
1M Tris-HCl	200	20mM
100% Triton X100	5	0.05%
Bring in solution at 40-60°C in water bath, takes about 30min then cool down on ice. Then add:		
Formamide	5500	55%
10% Blocking	1000	1%
Sigma Water	1500	

Buffer with 20% Formamide:

for CREN and EUR

Stock reagent	Volume (µl)	Final conc.
Dextran Sulfate	1g	10%
5M NaCl	1800	900mM
1M Tris-HCl	200	20mM
100% Triton X100	5	0.05%
Bring in solution at 40-60°C in water bath, takes about 30min then cool down on ice. Then add:		
Formamide	2000	20%
10% Blocking	1000	1%
Sigma Water	5000	

Amplification Buffer - 20ml:

Store at +4°C

Stock reagent	Volume (µl)	Final conc.
Dextran Sulfate	2g	10%
5M NaCl	8000	2M
10% Blocking	200	0.1%
1x PBS	11800	

DAPI mix - 2ml

Stock reagent	Volume (µl)	Final
DAPI 50µg/ml	40	1 µg/ml
1xPBS	140	0.5
Vectashield	280	1
Cititfluor	1540	5.5

DAPI mix - 5ml

Stock reagent	Volume (µl)	Final
DAPI 50µg/ml	200	2 µg/ml
1xPBS	350	0.5
Vectashield	700	1
Cititfluor	3750	5.5

Blocking Reagent 10% - 100ml

Prepare Maleic acid buffer

- 100mM Maleic acid
MW: 116.08g/mol
11.608g/l
- 150mM NaCl
MW: 58.44g/mol
8.766g/l
- adjust pH with a lot of solid NaOH to 7.5
- Add Boehringer Mannheim Blocking reagent to a final concentration of 10%
10g with 100ml Maleic Acid buffer
- Dissolve reagent on heating plate at around 60°C with stirring for about 1h, do not boil
- autoclave solution
- Aliquot to 10ml parts and store them frozen at -20°C; According to the manufacture solution should be good for a few days at 4°C and many month if stored frozen.

Paraformaldehyde 20%

- Turn on a 60 °C water bath in advance.
- Pour 150 mL of milliQ water into a 250 mL screw cap polycarbonate bottle. Add 50 g of paraformaldehyde.
- Fasten the cap securely and shake well to disperse the paraformaldehyde powder.
- Add 3 mL of 1 N, NaOH.
- Warm to 60 °C for 20 to 25 minutes shaking the bottle ever 5 minutes. Most but not all of the paraformaldehyde will dissolve in this time. It will still be quite cloudy. Don't worry. The last step is 0.2 µm filtration that will clean it up nicely.
- Add 3 mL of 1 N HCl
- Add 25 mL of 100 mM sodium phosphate buffer, pH 7.2.
- Add 25 mL of 5 M NaCl.
- Adjust pH to 7 using pH paper. This will take approximately 400 µL of 1N, NaOH. If you cook too long at 60°C it will be very acidic making it difficult to adjust the pH.
- Adjust the final volume to 250 mL with Milli Q using the graduations on the polycarbonate bottle.
- Filter using a 0.22 µm Millipore type GS filter (47 mm) using a designated formaldehyde side arm flask, frit and cup. There may be a substantial layer of paraformaldehyde on the filter.
- Transfer the filtrate to a clean 250 mL polycarbonate bottle. Label with the date concentration (20% formaldehyde) and your initials.
- Store at 4°C. Fresh formaldehydhe is good for 2 or days if kept in the refrigerator.

Phosphate buffers for PBS

10x PBS pH 7.7

2.76g NaH₂PO₄
14.24g Na₂HPO₄
75.97g NaCl
fill up to 1L with MQ

EDTA DisoDihy 0.5M

MW: 372.24g/mol
186.12g/l
add NaOH pastilles, only at pH 8 EDTA will dissolve!
fill up to 1 liter with MQ
pH 8.0; autoclave

SDS 10%

50g to 500ml MQ
do not autoclave!
dissolves at 40°C in water bath

Tris-HCL 1M

MW: 121.14g/mol
157.6g/l
With 6M HCL
pH 8.0
up to 1 Liter MQ

NaCL 5M

MW: 58.44g/mol
146.1g/500ml
autoclave

Abbreviations used

DMF:	N,N-Dimethylformamide $F_w=73.09\text{g/mol}$
TEA:	Triethylamine $F_w=101.19\text{g/mol}$
TYR-HCL:	Tyramine-HCL $M_w = 173.64\text{g/mol}$

Solutions

- DMF-TEA stock:**
1ml DMF + 10 μ l TEA
Prepare in 2ml Eppi
- TYR Stock:**
10mg TYR-HCL + 1ml DMF-TEA stock
100 μ l TYR.stock = 5.76 μ mol/l
Prepare in 15ml Greynner

Triethylamine & Dimethylformamide are dangerous: work in the hood and on ice!

- Succinimidyl ester:**
1mg active ester + 100 μ l DMF
1 mg Alexa₄₈₈ = 1.6 μ mol/l

Alexa ₄₈₈ :	$M_w = 643.41\text{g/mol}$ $\lambda_{\text{max}} = 495\text{nm}$ $E_{\text{max}} = 519\text{nm}$ $\epsilon = 71,000$ $C_f = 0.11$
------------------------	---

Alexa ₅₅₅ :	$M_w = 1250\text{g/mol}$ $\lambda_{\text{max}} = 555\text{nm}$ $E_{\text{max}} = 565\text{nm}$ $\epsilon = 150,000$ $CF_{280} = 0.08$
------------------------	---

Esters are light sensitive and prone to hydrolysis, therefore prepare shortly before tyramide synthesis; cool esters on ice until they are used for synthesis.

Synthesis

- 100 μ l Alexa₄₈₈ + 25.2 μ l TYR stock
100 μ l Alexa₅₅₅ + 12.6 μ l TYR stock
- Incubate at room temperature in the dark for ~12 hours, slowly rotating
- Dilute to 1ml with absolute Ethanol (874.8 μ l for A₄₈₈ and 887.4 μ l for A₅₅₅)
- Make aliquots (~50 μ l) and store at -20°C.

Chemicals needed

1mg A₄₈₈ Succinimidyl ester:
Molecular Probes Cat.Nr. A-20000

1mg A₅₅₅ Succinimidyl ester
Molecular Probes Cat.Nr. A-20009

N,N-Dimethylformamide (waterfree)
Sigma D 4551

Tyramine-HCL
Sigma T2879

Triethylamine
Sigma 17924

Microautoradiography

Cell sampling and fixing

1. Check the abundance of bacteria in your samples, to estimate the volume to filter (see step 6) and the amount of substrate to add. For example:

Abundance of bacteria in sample	Volume to filter (ml)
1×10^4 cell ml ⁻¹	50 (minimum)
5×10^5 cell ml ⁻¹	10
1×10^6 cell ml ⁻¹	5

2. Transfer samples in greyner tubes. 2 samples and 1 control. Kill control with PFA (final conc. ~2%) and wait for 15min.
3. Add radioactive labeled compounds, e.g.:

Substrate	Final concentration in sample*
³ H-Glucose	20 nM
³ H-Leucine	20 nM
³ H-Thymidine	20 nM
¹⁴ C-HCO ₃ ⁻	100 μCi

*Above concentrations are for marine systems

4. Incubate for at least 4h in the dark under in situ temperature. Take into account the expected level of activity for your sample and the substrate activity (e.g. for deep sea 24h for ³H-leucine, 72h for ¹⁴C-bicarbonate).
5. Kill all samples with PFA (final conc. ~2%) and incubate samples at 4°C in the dark for about 18h.
6. Filter the samples on 25mm 0.2 μm GTTP filters using a Manyfold and nitrocellulose support filters
7. Wash the filters 2x with 5ml fresh MQ.

Transfer of filter sections

All the following steps must be performed in the dark with a safety lamp, although complete darkness ensures a better result!

1. Melt photographic emulsion for 15min at 43°C in a water bath in the dark.
2. Coat the slides with the emulsion (see preparation notes) by dipping them into the tube (2 slide mailer).
3. Wipe away the emulsion from the underside of the slide with clean tissue paper.
4. Dry the slides for 5-10 min. (the time it takes to process 4 slides) on ice cold aluminum plate. This is best done in a small tray containing ice covered with aluminum foil.
5. One by one, take the prepared filter sections which are in nice order and carefully mount them on the slides, cells facing to the emulsion.
6. Once each slide is mounted, return it to the ice-tray, to be briefly dry, while processing the others
7. Place the slides in light proof box wrapped with aluminum foil, containing silica gel drying beads, for exposure at 2- 4°C. The exposure time is dependent on the activity of the substrate and expected activity of the cells. See appendix for table of suggested guidelines.

*For substrates not included in the table an appropriate exposure time should be tested.

Developing and Fixing

1. Place box with slides in the dark room (17°C) for one hour before starting.
2. Prepare developing chemicals in slide staining jars

All the following steps must be performed in complete darkness!

3. Remove aluminium foil and place the slides in the slide rack.
4. Develop the slides using Kodak specifications: 2 min in Dektol Developer, 10s to stop developing in MQ, 5 min in fixer, 2 min in MQ.

*The Dark room should be at around 17°C, especially the developer and fixer. All developing chemicals and MQ should be kept at this room temperature.

After this you can turn on the light!

1. Keep the slides in the MQ until you are back in the lab. Dry the underside of the slide
2. Surround the borders of the filters on the underside of the slide with a marker. This is also the last chance to know from the filter labels which section is which - so label the slide accordingly.
3. Dry a bit with paper the side with filters, but not too much, otherwise it is difficult to remove the filters. If the filter is too dry and sticks to the slide, it can be rewet with a drop of water.
4. Carefully peel the filter sections off
5. Put a drop of DAPI-mix in each filter outline and put a cover slip over it.
6. They can be stored at -20°C until they are counted on the epifluorescence microscope.

Preparing the emulsion:

The Emulsion is a mix of NTB and Sigma water 2:1 (v/v) = 7ml + 3.5ml

1. Melt the emulsion for 45-60min at 43-45°C before use.
2. Wrap the 2 slide mailer (tubes) with colored PVC tape and ensure that the lid is covered with black marker pen.
3. Add 3.5 ml Sigma water.
4. In the dark room, add 7ml of photographic emulsion.

The emulsion is highly sensitive to light, therefore, complete darkness is recommended. This is best done using a bottle top dispenser (e.g. Dispensette, Sigma-Aldrich).

5. Wrap the tubes in aluminum foil and store at 4°C.

*Take into account the emulsion has an expiration date of few months. Each tube with the emulsion can be used a maximum of three times.

Suggested exposure time for various substrates according to literature

Substrate	Concentration (nM)	Emulsion	Incubation Time (h)	Temperature (°C)	Location	Reference
³ H-AA mixture (47 Ci/mmol ⁻¹)	2.1	NTB-2	48	-20	Delaware Bay	Cottrell and Kirchman 2000
³ H-NAG (9,9 Ci/mmol ⁻¹)	10	NTB-2	48	-20		
³ H-Protein [from <i>Vibrio</i> grown on ³ H-leucine]		NTB-2	48	-20		
³ H-Chitin [from fungus grown on ³ H-NAG]		NTB-2	48	-20		
³ H-Thymidine (83,5 Ci mmol ⁻¹)	20	NTB-2	48	4	Delaware estuary	Cottrell and Kirchman 2003
³ H-Leucine (150 Ci mmol ⁻¹)	20	NTB-2	48	4		
³ H-AA mixture (47 Ci mmol ⁻¹)	0.5	NTB-2	48-144	4	Noth Atlantic Ocean	Malmstrom et al. 2004
³⁵ S-DMSP (1,170 Ci mmol ⁻¹)	<0.1	NTB-2	336-480	4		
³ H-Thymidine (83,5 Ci mmol ⁻¹)	20	NTB-2	28-48	4	Delaware estuary	Cottrell and Kirchman 2004
³ H-Leucine (150 Ci mmol ⁻¹)	20	NTB-2	3-7	room T		
³⁵ S-DMSP (12-43 TBq mmol ⁻¹)	<0.1	LM-1	48-144	4	Noth Atlantic Ocean	Malmstrom et al. 2004
³ H-Protein (0,3-0,9 Ci/gr ⁻¹) [from <i>Vibrio</i> grown on ³ H-Leucine]	20 ng/ml	NTB-2	10 days	4	Noth Atlantic Ocean	Malmstrom et al. 2005
³ H-AA mixture (47 Ci/mmol ⁻¹)	0.5	NTB-2	48-72	4		
³ H-Glucose (40 Ci/mmol ⁻¹)	0.5	NTB-2	48-72	4		
³ H-Glucose (33 Ci/mmol ⁻¹)	2	NTB-2	12-24	4	Delaware estuary	Elifantz et al. 2005
³ H-EPS [from <i>Nitzschia</i> grown on ³ H-glucose]	1,5	NTB-2	72-144	4		
³⁵ S-DMSP (12-43 TBq/mmol ⁻¹)	<0.1	LM-1	48-144	4	North Atlantic Ocean	Malmstrom et al. 2005
³⁵ S-MeSH (43 TBq/mmol ⁻¹)	<0.01	LM-1	48-144	4		

NAG: N-acetylglucosamine

EPS: extracellular polymeric substances

DMSP: dimethylsulfoniopropionate

MeSH: methanethiol

Product	Description	Company	Art.Nr.:	Size
<u>Agarose</u>	low melting point sea plaque GTGagarose	FMC bioProducts	56111	
<u>Aspartic acid</u>	D-[2,3 ³ -H]	Amersham Bioscience	Nr.: TRK606	250µCi
<u>Aspartic acid</u>	L-[2,3 ³ -H]	Amersham Bioscience	Nr.: TRK574	1mCi
<u>Blocking Reagent</u>	Boehringer Mannheim Blocking reagent	Roche Diagnostics GmbH (NL)	1-096-176	50g
<u>Buechner funnel</u>	(diam.59mm)	VWR International	HALD127C/1	
<u>Citifluor</u>	Glycerol-PBS solution	Citifluor Ltd.	AF1	25ml
<u>DAPI</u>	DAPI stain	Sigma-Aldrich Chemie BV	D9564	1x10mg
<u>Dextran Sulfate</u>		Sigma-Aldrich Chemie BV	D8906	100g
<u>Dispensette</u>	0.5-5ml	Sigma-Aldrich Chemie BV (NIOZ:Brand)	Z33,194-5 (07 X5250?)	
<u>DMF</u>	NN-Dimethyl formamide	Sigma-Aldrich Chemie BV	D4551	250ml
<u>EDTA</u>		Plusone	17-1324-01	
<u>FISHprobes</u>		Thermo Electron Corporation Biopolymers		
<u>Formamide</u>	(Fluka)	Sigma-Aldrich Chemie BV	47671	1l 250ml
<u>Glascrobe</u>		Bel-Art Products, Inc. Scienceware	H441500000	
<u>H₂O₂</u>	30% (w/w)	Sigma-Aldrich Chemie BV	H1009	5ml
<u>Kodak Dektol Developer</u>		Sigma-Aldrich Chemie BV	P6917	1GA
<u>Kodak Fixer</u>		Sigma-Aldrich Chemie BV	P6557	1GA
<u>Kodak NBT-2</u>		Integra Biosciences GmbH	IB-01433	1x
<u>Lysozyme</u>	from chicken egg,white	Sigma-Aldrich Chemie BV	Nr.: L7651	50000units/mg 25g
<u>maleic acid</u>		Sigma-Aldrich Chemie BV	M-0375	500g
<u>multiwell plates</u>	12-well (TC treated with lid)	Sigma-Aldrich Chemie BV	M 8687	1St.
<u>Nitril-gloves</u>	Format (blue)	Unigloves- Firma Kleiss(in AUT)	Art.Nr.: 6003(M)	100St

Product	Description	Company	Art.Nr.:	Size
<u>Nitrocellulosefilter</u>	0.45µm25mm(diam.)	Sigma-AldrichChemie BV	N9020	100EA
<u>Paraformaldehyde</u>		Sigma-Aldrich Chemie BV	Nr.: P6148	1x500g
<u>Polycarbonate filter</u>	0.2µm 25mm(diam.)	Sigma-Aldrich Chemie BV	P9199	100
<u>Proteinase K</u>	(Fluka) from triturachium album	Sigma-Aldrich Chemie BV	Nr.: 82456	1x1ml
<u>Safetylight</u>	Dunkel kammer leuchte 230V,50Hz	Kaiser Fototechnik	Nr.: 4018	
<u>Safetylight</u>	lamp model B	Eastman Kodak Company	141-2212	
<u>Safetylight filter</u>	Dunkel kammer filter (rot) 9x12cm	Hama	Nr.:8194	
<u>Safetylight filter</u>	Kodak 2 darkred	Eastman Kodak Company	152-1525	
<u>SDS</u>	Lauryl Sulfate	Sigma-Aldrich Chemie BV	L-4509	100g
<u>2 slide mailer</u>	tubes for emulsion	Raymond A Lamb Limited	Nr.:E6.2	125x
<u>staining dishes</u>		Sigma-Aldrich Chemie BV	S4642	2x3
<u>staining dishes</u>	slide rack	Sigma-Aldrich Chemie BV	S5017	1x3
<u>staining dishes</u>	rack handle	Sigma-Aldrich Chemie BV	S5142	1x3
<u>TEA</u>	Triethyl amine	Sigma-Aldrich Chemie BV	17924	
<u>TRIS</u>		Plusone	17-1321-01	
<u>Triton X100</u>	X-100	Sigma-Aldrich Chemie BV		500ml
<u>Tyramide dye</u>	Alexa Fluor 488 succinimidyl ester	Molecular Probes Europe BV	Nr.: A-20000	1x1mg
<u>Tyramide dye</u>	Alexa Fluor 555 succinimidyl ester	Molecular Probes Europe BV	Nr.: A-20009	1x1mg
<u>TYR-HCl</u>	Tyramine- HCl	Sigma-Aldrich Chemie BV	T2879	5g
<u>Vectashield</u>		Vector Laboratories, Inc.	H-1000	10ml
<u>Water</u>	Sigma Water	Sigma-Aldrich Chemie BV	W4502	1L